

# Package ‘MetabolomiQCsR’

June 9, 2017

**Type** Package

**Title** QC for LC-MS data

**Version** 1.3

**Date** 2017-06-09

**Author** Jan Stanstrup

**Maintainer** Jan Stanstrup <stanstrup@gmail.com>

**Description** Functions for performing QC on LC-MS.

**Depends** R (>= 3.0)

**Imports** tibble, magrittr, dplyr, purrr, xcms, messageR, tidyr, scales, gg-plot2, ini, stringr, RCurl, readr, utils, WGCNA, pool, RMySQL, DBI, lubridate, httr, purrrlyr

**Suggests** knitr, rmarkdown, faahKO, plotly, svglite, chemhelper, RColorBrewer, heatmaply, listviewer

**License** GPL (>= 2)

**URL** None

**biocViews** MassSpectrometry, Metabolomics

**VignetteBuilder** knitr

**RoxygenNote** 6.0.1

**Remotes** stanstrup/messageR, stanstrup/chemhelper, rstudio/pool

## R topics documented:

closest_match . . . . .	2
EIC_contaminants . . . . .	2
extract_polarity . . . . .	3
get_cont_list . . . . .	4
get_EICs . . . . .	4
ic_boxes . . . . .	5
ic_measurements . . . . .	6
ic_token . . . . .	6
parse_filenames . . . . .	7
peak_factor . . . . .	8
plotly_clean_tt . . . . .	8
plot_chrom . . . . .	9
plot_contaminants . . . . .	9
tbl2ROI . . . . .	10
xcmsRaw_to_tbl . . . . .	10

**Index****11**


---

closest_match	<i>Match list of standards to peak table</i>
---------------	--

---

**Description**

This function will match a table of standard compounds and a peak table by m/z and retention time. If there is more than one possible hit the highest intensity peak will be chosen.

**Usage**

```
closest_match(stds, peakTable, rt_tol = 0.25, mz_ppm = 30, rt_col = "rt",
              mz_col = "mz", int_col = "into")
```

**Arguments**

stds	<a href="#">tibble</a> of standards to match to a peak table
peakTable	<a href="#">tibble</a> containing peak table supplied by <a href="#">findPeaks</a> (but converted to <a href="#">tibble/data.frame</a> ).
rt_tol	Retention time tolerance for matching peaks. Pay attention to the unit of your tables. rt_tol should match and stds and peakTable should use same units (i.e. minutes or seconds).
mz_ppm	ppm for matching peaks.
rt_col	Character string giving the column containing the retention times. Must be same in standards and peak table.
mz_col	Character string giving the column containing the m/z values. Must be same in standards and peak table.
int_col	Character string giving the column containing the intensities in the peak table.

**Value**

A vector having the length equivalent to the number of rows in stds giving the indices of the hits in peakTable.

---

EIC_contaminants	<i>Find EICs that behave like contaminants</i>
------------------	--

---

**Description**

This function looks in raw LC-MS data for "features"/EICs that behave like contaminants. Behaving like contaminants in this case means that a certain m/z values is present in more than min\_time above intensity min\_int.

**Usage**

```
EIC_contaminants(raw, bin_ppm = 30, interval_ppm = 30, min_time = 5,
                 merge_corr = 0.9, merge_ppm = 30, min_int = 5000)
```

### Arguments

raw	<a href="#">xcmsRaw</a> object to profile
bin_ppm	Tolerance (ppm) for initial binning of m/z values
interval_ppm	Tolerance for creating final EICs after merging similar bins
min_time	Minimum time (minutes) an EIC should be above min_int to be considered a contaminant.
merge_corr	Minimum correlation between EICs to be merged.
merge_ppm	Maximum difference (ppm) between EICs to be merged.
min_int	Minimum intensity that the EIC needs to be above for a minimum of min_time.

### Value

A [tibble](#) containing the columns:

- **mz:** m/z of the proposed contaminant
- **EIC:** EIC of the m/z.

---

extract_polarity	<i>Extract polarity from xcmsRaw object.</i>
------------------	--

---

### Description

Extracts polarity from an xcmsRaw object. The polarity found in the majority of scans is returned.

### Usage

```
extract_polarity(xraw)
```

### Arguments

xraw	The xcmsRaw object to extract polarity from.
------	--

### Value

A character string giving the polarity. Can be "positive", "negative", or "unknown".

---

get_cont_list	<i>Get list of known contaminants</i>
---------------	---------------------------------------

---

### Description

Get list of known contaminants

### Usage

```
get_cont_list(polarity = c("positive", "negative", "unknown"), type = "URL")
```

### Arguments

polarity	The polarity to get contaminants for. Can be "positive", "negative" or "unknown". If "unknown" the list specified in the MetabolomiQCs.conf is used. MetabolomiQCs.conf can be in the working folder or the home folder. If those are not found the package default is used (unknown will used the positive mode list).
type	If using local or remote. Only "URL" implemented which downloads a list from <a href="https://github.com/stanstrup/common_mz">https://github.com/stanstrup/common_mz</a>

### Value

tbl A [tibble](#) containing the columns:

- **Monoisotopic ion mass (singly charged):** m/z of the contaminant
- **Ion type:** Notation for adduct/fragment type
- **Formula for M or subunit or sequence:** Molecular formula
- **Compound ID or species:** Name of the compound
- **Possible origin and other comments:** Suggestion for the origin of the contaminant
- **References:** Reference for the contaminant

### References

[https://github.com/stanstrup/common\\_mz](https://github.com/stanstrup/common_mz)

---

get_EICs	<i>Get EICs from xcmsRaw object</i>
----------	-------------------------------------

---

### Description

Takes an [xcmsRaw](#) object and extracts EICs. Can do multiple ranges and exclude certain masses unlike [getEIC](#). Can be used to extract the TIC too.

### Usage

```
get_EICs(xraw, range_tbl, exclude_mz = NULL, exclude_ppm = 30,
  range_tbl_cols = c("mz_lower", "mz_upper"), BPI = FALSE, min_int = 0)
```

### Arguments

xraw	<code>xcmsRaw</code> object to get EIC(s)/TIC from.
range_tbl	data.frame/ <code>tibble</code> with columns for the lower and upper m/z boundaries of EIC slice(s).
exclude_mz	Masses to exclude from the EIC. Most useful to remove contaminants from TICs.
exclude_ppm	ppm tolerance of exclude_mz
range_tbl_cols	Which columns in range_tbl holds the lower and upper range. defaults to c("mz_lower", "mz_upper").
BPI	Logical selecting to calculate TIC (FALSE) or BPI.
min_int	the minimum intensity mass peak to include

### Value

tbl A `tibble` containing the columns:

- **scan:** scan number
- **scan\_rt:** Retention time of scan
- **intensity:** The summed intensity for each scan in the given m/z interval

---

ic_boxes	<i>Convert raw data into a tibble of xcmsRaw objects.</i>
----------	---

---

### Description

Convert raw data into a tibble of `xcmsRaw` objects.

### Usage

```
ic_boxes(token)
```

### Arguments

token	A character vector containing the ICMeter token
-------	---

### Value

data.frame containing information about each box registered on the token account

---

ic_measurements	<i>Convert raw data into a tibble of xcmsRaw objects.</i>
-----------------	---

---

**Description**

Convert raw data into a tibble of xcmsRaw objects.

**Usage**

```
ic_measurements(token, boxQR, fromDate, toDate)
```

**Arguments**

token	A character vector containing the ICMeter token
boxQR	A character vector containing QR code for the box you want to query
fromDate	The starting date to draw data from
toDate	The ending date to draw data to. A maximum of one month can be draw at a time.

**Value**

data.frame holding all the data drawn

---

ic_token	<i>Convert raw data into a tibble of xcmsRaw objects.</i>
----------	---

---

**Description**

Convert raw data into a tibble of xcmsRaw objects.

**Usage**

```
ic_token(user, pass)
```

**Arguments**

user	ICMeter username
pass	ICMeter password

**Value**

A character vector containing the ICMeter token

---

parse_filenames	<i>Parse filenames from given string</i>
-----------------	--

---

## Description

Parse filenames from given string

## Usage

```
parse_filenames(filenames, filenamemask)
```

## Arguments

filenames	Character string with the filenames to parse.
filenamemask	Filename mask. Variables to capture are inclosed by "%". See examples

## Value

data.frame with fields parsed from the filenames

## Author(s)

Stinus Lindgreen, <sldg@steno.dk>; Jan Stanstrup, <jpzs@steno.dk>

## Examples

```
testfilenames <- c("RemoveThis_METSY_sample1___17-random_AndThis.txt",
  "RemoveThis_METSY_sample2___21-something_AndThis.txt",
  "RemoveThis_METSY_sample3___123-blah_AndThis.txt",
  "RemoveThis_JDRF_SameName___1-dummyA_AndThis.txt",
  "RemoveThis_JDRF_SameName___2-dummyB_AndThis.txt",
  "RemoveThis_JDRF_SameName___3-dummyC_AndThis.txt",
  "RemoveThis_this_file___name-is_wrong",
  "RemoveThis_and_this_AndThis.txt",
  "RemoveThis_this_one___should-work",
  "RemoveThis_how_about___this_-_one_AndThis.txt"
)

defaultmask <- "RemoveThis_%study%_%name%___%rep%-%dummy%_AndThis.txt"

parse_filenames(testfilenames,defaultmask)
```

---

peak_factor	<i>Calculate Tailing Factor and Asymmetry Factor</i>
-------------	--

---

**Description**

Calculate Tailing Factor and Asymmetry Factor

**Usage**

```
peak_factor(EIC, rt, factor = "TF")
```

**Arguments**

EIC	EIC containing the peak to calculate for. <a href="#">tibble</a> as produced with <a href="#">get_EICs</a> .
rt	Retention time of the center of the peak (Numeric)
factor	to calculate. Character string either "TF" (Tailing Factor) or "ASF" (Asymmetry Factor).

**Value**

Numeric

**References**

<http://www.chromforum.org/viewtopic.php?t=20079>

---

plotly_clean_tt	<i>Make replacements in plotly tooltips</i>
-----------------	---

---

**Description**

Make replacements in plotly tooltips

**Usage**

```
plotly_clean_tt(plotly, rep)
```

**Arguments**

plotly	A <a href="#">plotly</a> object.
rep	A named character vector. Names are the text to replace and the string is the replacement string.

**Value**

A [plotly](#) object.



---

plot_chrom	<i>Plot chromatogram</i>
------------	--------------------------

---

**Description**

Plot chromatogram

**Usage**

```
plot_chrom(tbl, RT_col = "RT", Intensity_col = "Intensity")
```

**Arguments**

tbl	tbl with retention time and intensity to plot
RT_col	Name of the retention time column
Intensity_col	Name of the intensity column

**Value**

a [ggplot](#) object.

---

plot_contaminants	<i>Bar plot of contaminants</i>
-------------------	---------------------------------

---

**Description**

Bar plot of contaminants

**Usage**

```
plot_contaminants(data, title, x_var = "comp_name", y_var = "EIC_median")
```

**Arguments**

data	tbl with the contamination amounts
title	Plot title
x_var	Column name that holds the compound/contaminant names
y_var	Column name that holds the compound/contaminant values

**Value**

a [ggplot](#) object.

---

tbl2ROI	<i>Convert a list of peaks (rt / m/z pairs) to a Region of Interest (ROI) list for use with <a href="#">findPeaks</a></i>
---------	---

---

### Description

Convert a list of peaks (rt / m/z pairs) to a Region of Interest (ROI) list for use with [findPeaks](#)

### Usage

```
tbl2ROI(tbl, raw, ppm, rt_tol)
```

### Arguments

tbl	<a href="#">tibble</a> containing the columns "rt" and "mz". rt needs to be in seconds.
raw	<a href="#">xcmsRaw</a> object to create ROI for. It needs to be a specific <a href="#">xcmsRaw</a> to match retention times to scan numbers.
ppm	ppm tolerance for the generated ROI.
rt_tol	Retention time tolerance (in sec!) for the generated ROI.

### Value

List containing the ROIs. Each list contains mz, mzmin, mzmax, scmin, scmax, length (set to -1, not used by centWave) and intensity (set to -1, not used by centWave) columns.

---

xcmsRaw_to_tbl	<i>Convert raw data into a tibble of xcmsRaw objects.</i>
----------------	---

---

### Description

Convert raw data into a tibble of xcmsRaw objects.

### Usage

```
xcmsRaw_to_tbl(files, ...)
```

### Arguments

files	character vector of file names/paths.
...	further arguments to <a href="#">xcmsRaw</a> .

### Value

A [tibble](#) containing the columns:

- **file:** Filename without path.
- **polarity:** Character string of "positive", "negative", or "unknown".
- **raw:** The xcmsRaw objects.
- **path:** The input path (files).

# Index

`closest_match`, [2](#)

`data.frame`, [2](#)

`EIC_contaminants`, [2](#)  
`extract_polarity`, [3](#)

`findPeaks`, [2](#), [10](#)

`get_cont_list`, [4](#)  
`get_EICs`, [4](#), [8](#)  
`getEIC`, [4](#)  
`ggplot`, [9](#)

`ic_boxes`, [5](#)  
`ic_measurements`, [6](#)  
`ic_token`, [6](#)

`parse_filenames`, [7](#)  
`peak_factor`, [8](#)  
`plot_chrom`, [9](#)  
`plot_contaminants`, [9](#)  
`plotly`, [8](#)  
`plotly_clean_tt`, [8](#)

`tbl2ROI`, [10](#)  
`tibble`, [2–5](#), [8](#), [10](#)

`xcmsRaw`, [3–5](#), [10](#)  
`xcmsRaw_to_tbl`, [10](#)